

Helpful suggestions and top considerations for good qPCR

- 1.) RNA quality/integrity assessments, e.g. achieving an “RNA Integrity Number” (RIN) of 7 and above for all RNA samples (e.g. using an Agilent Bioanalyzer 2100) is an absolute necessity preceding downstream activities including reverse transcription and subsequent qPCR using the resulting cDNA. See <http://dna.biotech.iastate.edu/nextgen.html> for more info on required RNA and DNA quality and quantity assessments.
- 2.) Splice variants exist for >60% of all known mRNA transcripts in most species, so, during primer-probe design, it is of utmost importance to acknowledge this. Know which splice variant(s) you are targeting with each of your primer-probe designs (or target them all at the same time). Check all designs against all available data bases including BLAST (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>) and BLAT (<http://genome.ucsc.edu/cgi-bin/hgBlat>). In the event that splice variants are a non-issue, pseudo-genes and other homologous sequences must be screened for as well in order to be certain that chosen primer-probe(s) are not inadvertently sequestered away from the intended final target PCR/qPCR reactions - rendering them unacceptably less efficient.
- 3.) Choice of Reverse Transcriptase enzyme for cDNA synthesis is a major decision as some RT enzymes are very poor (e.g. Multiscribe™ RT) at yielding plentiful and full-length cDNAs. Suggested RTases are: Quanta BioScience’ QuantaScript™, Invitrogen’s SuperScript III™ and Takara’s PrimeScript™ for first strand synthesis (e.g. for Two-Step qPCR).
- 4.) Ample dilution of RNA samples beyond their RT-inhibitory range, and dilution from there, to the same ng/μl preceding their addition to RT reactions is an essential practice for generating cDNA intended for subsequent use in qPCR.
- 5.) Priming of RT reactions. A mixture of oligo d[T] and random hexamers, octamers or nonamers has shown to be the most effective, although sole use of random pentadecamers has been lauded as even more effective by some. Oligo d[T] priming is not suggested for RNA extracted from formalin-fixed tissue, as many RNAs lose their poly A tails after such tissue treatment. If oligo d[T] priming is used, 3’-bias is a concern, therefore, primer-probes should be designed to target sequences nearer the 3’ end rather than the 5’ end of the known mRNA transcripts. 18S rRNA and histones mRNAs do not contain poly A tails from the start, so oligo d[T] priming is not appropriate for those targets – although mis-priming by some RT enzymes can yield varying amounts of 18S rRNA cDNA when reverse transcription is carried out using oligo d[T] for first-strand priming.
- 6.) Using the “PREXCEL-Q Method” is suggested for all qPCR approaches. All dilution-related complications of qPCR are dealt with using this tool. And this tool has been endorsed by some of the world’s foremost qPCR experts (Bustin, Vandesompele).
- 7.) Removal of inhibitors from DNA samples (e.g. using MO Bio Labs reagents) is highly suggested.